

# Digital refractometers DR7000 series, with audit trail function

-Compliant with CFR 21 Part 11

|                     |                                  |
|---------------------|----------------------------------|
| Dimensions (WxDxH): | 230 x 360 x 170 mm               |
| Weight:             | 5 kg                             |
| Power supply:       | 90 ... 240 V&middot;AC, 50/60 Hz |
| IP code:            | IP 20                            |

| Type     | Measured range                       | Resolution           | Accuracy                             | Min. temperature<br>°C |
|----------|--------------------------------------|----------------------|--------------------------------------|------------------------|
| DR7200-P | nD 1.3000 ... 1.58000 ... 95 %Brix   | nD 0.00010.1 %Brix   | nD &plusmn;0.0001&plusmn;0.1 %Brix   | 10                     |
| DR7400-P | nD 1.30000 ... 1.720000 ... 95 %Brix | nD 0.000010.01 %Brix | nD &plusmn;0.00002&plusmn;0.02 %Brix | 10                     |

Zubehör:  
Accessories for digital refractometers DR7000 series / A. Kr]ss Optronic