

Clear Seal 3730

Thin clear Polyester backed film forming a piercable seal to PP, PS and COC plates.

- Application: ABi 3730 capillary sequencing. Quantitative PCR. Not peelable.
- Sealing: 1 sec, 170 to 175 °C
- Seal Integrity: -80 to 80 °C, 110 °C with PCR lid pressure
- 1000 sheets on perforated roll

Sterile: No

Material	Colour	Pierceable	Self-adhesive	Sealable	PK	Cat. No.
PE	Transparent	Yes	No	Yes	1	9.776 009