

Semi micro cuvettes

Capacity: 1.5 ... 3 ml
Ext. dim. (WxDxH): 12.5 x 12.5 x 45 mm
Type of cuvette: Semi-Micro

Wavelength nm	Path length mm	Min. capacity ml	Max. capacity ml	Material	Ext. width mm	Ext. depth mm	Ext. height mm	PK	Cat. No.
300 ... 900	10	1.5	3	PMMA	12.5	12.5	45	100	9.406 115
340 ... 900	10	1.5	3	PS	12.5	12.5	45	100	9.406 114