

Vacuum measuring instrument DCP 3000, without sensor

Dimensions (WxDxH): 138 x 124 x 114 mm
Weight: 440 g
Plug type: EU/UK/CH/US/AUS
Power supply: 100 ... 230 V, 50/60 Hz

Type	Width mm	Depth mm	Height mm	Type of auxiliary energy	PK	Cat. No.
DCP 3000	138	124	114	Mains connection	1	7.940 260

Zubehör:
Accessories for Vacuum Controller VACUU-SELECT and Measuring instrument DCP 3000 / Vacuubrand
Accessories for Vacuum Controller VACUU-SELECT and Measuring instrument DCP 3000 / Vacuubrand