

Chromatographic columns, borosilicate glass

NS necks are in compliance with DIN 12242 standards.
Material: Borosilicate
With stopcock: Yes
Material stopcock: PTFE key
Standard: DIN 12242

Ø mm	Length mm	Dimensions (ØxL)	Socket[s]	PK	Cat. No.
10	200	10 x 200 mm	NS14/23	1	4.008 398
20	400	20 x 400 mm	NS29/32	1	4.008 399
30	600	30 x 600 mm	NS29/32	1	4.008 400

Alternativen:
Chromatographic columns, DURAN® tubing / Lenz-Laborglas