

Dual channel pH/ion meter Orion®; Pro Star PI217

Scope of supply: Dual channel pH/ion meter, electrode stand, USB flash drive,
PC cable, universal power supply with plug adapter EU, UK, US, AUS and China
pH range: -2.000 ... 20.000
pH range resolution: 0.1/0.01/0.001
pH range accuracy: ± 0.002
mV range: -2000 ... 2000
mV range accuracy: ± 0.2 mV or ± 0.05 % of the reading (whichever is greater)
Ion concentration (ISE): 0.0001 ... 19999
Ion concentration (ISE) resolution: 0.0001 minimum, 1 to 5 digits
Ion concentration (ISE) accuracy: ± 0.2 or ± 0.05 % of the reading (whichever is greater)
Temperature: -30 ... 130 °C
Temperature resolution: 0.1 °C
Temperature accuracy: ± 0.3 °C
Interfaces: 2 x USB
Dimensions (WxDxH): 153 x 219 x 92 mm
Weight: 1 kg
Power supply: 100 ... 240 V, 50/60 Hz
IP code: IP 54
Warranty: 2 years

Type	Min. pH value	Max. pH value	Width mm	Depth mm	Height mm	Type of auxiliary energy	PK	Cat. No.
Orion®; Pro Star PI217	-2	20	153	219	92	Mains connection	1	4.711 698