

Soxhlet extractors, without stopcock

Borosilicate glass 3.3 which is resistant to heat and almost all chemicals. They meet the highest quality standards according to DIN 12602.
Standard: DIN 12602

Nominal capacity ml	Cone[s]	Socket[s]	Material	With stopcock	PK	Cat. No.
30	NS29/32	NS29/32	Boro 3.3	No	1	4.008 417
70	NS29/32	NS34/35	Boro 3.3	No	1	4.008 418
100	NS29/32	NS45/40	Boro 3.3	No	1	4.008 419
150	NS29/32	NS45/40	Boro 3.3	No	1	4.008 420
250	NS29/32	NS45/40	Boro 3.3	No	1	4.008 421
500	NS29/32	NS45/40	Boro 3.3	No	1	6.243 673

- Zubehör:
PTFE sealing rings for jointed glassware / Dr. Dietmar Glindemann
Sleeves, cone joint, PTFE / Bohlender
Sleeves with ribs, PTFE / Bohlender
Joint clips KECK, POM, for conical ground joints / DWK Life Sciences
Joint clips KECK, POM, for spherical ground joints / DWK Life Sciences
Ground joint clips for sleeve connections, wire, Chrome-nickel steel / H W S Hans W. Schmidt
Joint clips for conical joint sleeve connections, stainless steel, with set screw / Lenz-Laborglas
Joint grease, K.W.S / J.P. PÅllath
BUDDE - joint grease, silicone-free / Buddeberg
Engler distilling flasks, DURAN®, ASTM, with beaded rim / DWK Life Sciences
Gas leak detection fluid, bottle / Otto Fahsig