

96-Well Optical-Bottom plates, with CC2™ surface

Achieve optimal results in cell-based fluorescence and luminescence applications using the Thermo Scientific™ Nunc™ 96-Well Optical-Bottom plates with coverglass base. Signals are intensified through the highly optically clear coverglass bottom due to improved signal-to-noise ratios.

Material:	PS
Colour:	Black
Bottom shape:	Flat
Growth surface:	0.32 cm²
Sterile:	Yes
With lid:	Yes

No. of wells	Nominal capacity µl	PK	Cat. No.
96	400	30	4.009 176