

# Magnetic stirrer MSL 50 Digital

Max. stirring capacity: 50 l  
Speed: 20 ... 1,000 1/min  
Material stirring surface: Stainless steel  
Material housing: Metal/Epoxy  
Permitted ambient temperature: 5 ... 40 °C  
Dimensions (WxDxH): 354 x 369 x 77 mm  
Weight: 4 kg  
Rated capacity: 10 W  
Power supply: 100 ... 240 V, 50/60 Hz

| Type           | No. of stirring positions | Min. speed 1/min | Max. speed 1/min | Display | Min. permitted ambient temperature °C | Max. permitted ambient temperature °C |
|----------------|---------------------------|------------------|------------------|---------|---------------------------------------|---------------------------------------|
| MSL 50 Digital | 1                         | 20               | 1000             | Digital | 5                                     | 40                                    |
| MSL 50 Digital | 1                         | 20               | 1000             | Digital | 5                                     | 40                                    |

Alternativen:  
Magnetic stirrer MSL 25 Digital / Velp Scientifica SRL