

Desiccators Maxi 1 Protect, PMMA

Aluminium frame. Orange, non-transparent acrylic glass panels reduce light incidence and protect against ultraviolet radiation. Including four shelves made of stainless steel, tray and hygrometer, four castors, usable with up to 17 shelves.
Int. Dim. (WxDxH): 495 x 540 x 1,030 mm
Ext. dim. (WxDxH): 560 x 580 x 1,150 mm
Max. load shelves: 30 kg
Weight: 30 kg
Temperature resistance: -20 ... 70 °C
Material: PMMA

Nominal capacity l	Max. no. of shelves	Int. width mm	Depth int. mm	Height int. mm	Ext. width mm	Ext. depth mm	Ext. height mm	Min. t. resista °C
311	34	495	540	1030	560	580	1150	-20