

Bottle carriers, wire/PE

Made of PE coated wire.
For carrying laboratory bottles of various sizes.
With base and fixed carrying handle.

For bottles ml	No. of positions	Compartment Ø mm	Height mm	Overall height mm	Colour	Material	PK	Cat. No.
250	4	80	80	320	White	Wire	1	9.102 224
250	6	80	80	320	White	Wire	1	9.102 226
250	8	80	80	320	White	Wire	1	9.102 228
250	10	80	80	320	White	Wire	1	9.102 230
250	12	80	80	320	White	Wire	1	9.102 232
1000	4	100	100	330	White	Wire	1	9.102 214
1000	6	100	100	330	White	Wire	1	9.102 216
1000	8	100	100	330	White	Wire	1	9.102 218
1000	10	100	100	330	White	Wire	1	9.102 220
2000	4	130	110	370	White	Wire	1	9.102 204
2000	6	130	110	370	White	Wire	1	9.102 206
2000	8	130	110	370	White	Wire	1	9.102 208

Alternativen:
Bottle carriers, wire/plastic coated / W]rtt. Drahtwarenfabrik
Bottle carriers for Duran square bottles / W]rtt. Drahtwarenfabrik