

Roller bottles WHEATON^{®}, borosilicate glass, with screw cap, PP

For adherent laboratory cell culture. Meets the requirements of USP Type I and ASTM E-438 Type I, Class A. Autoclavable. White screw cap 45 mm, with inner sealing ring and pouring ring.
Material:

Boro 3.3							
Nominal capacity ml	Ø mm	Height mm	Dimensions (ØxH)	Material lid	With lid	PK	Cat. No.
1380	110	240	110 x 240 mm	PP	Yes	4	4.662 207
1760	110	285	110 x 285 mm	PP	Yes	4	4.662 208