

384 Well LoBase microplates, Small Volume™, PS, with transparent bottom

LoBase microplates can be used for bottom reading measurements of samples with low working volume. The plates have a transparent µ;Clear[®] bottom. Due to the low level of autofluorescence, the plates are well suited for cell-based test systems. Free of human DNA, RNase/DNase, and pyrogens.

Working volume:	4 ... 25 µ;l
Material:	PS
Colour:	Black
Bottom shape:	Flat
Coating:	Non-treated
Sterile:	No
With lid:	No

No. of wells	Nominal capacity µl	Autoclavable	RNase/DNase-free	Pyrogen free	PK	Cat. No.
384	28	No	Yes	Yes	80	6.286 346

Alternativen:
96 Well Polystyrene Microplates, PS / Greiner Bio-One
384 well microplates, PS, non-treated / Thermo Elect.LED
384 Well HiBase microplates, Small Volume™, PS / Greiner Bio-One