

Desiccator Control-Maxi 2

Desiccator with aluminium frame, acrylic glass panes, and castors. Two separate chambers with 2 stainless steel shelves per chamber. The desiccator is stackable and can be equipped with up to 8 shelves per chamber. The automatic control unit continuously regulates the air humidity to the setpoint. The setpoint can be digitally adjusted from 15 to 95 % RH. No manual filling or emptying of water necessary. The display shows the temperature, the current operating status (dehumidifying/humidifying), as well as the set target value, and current actual value of the air humidity. Internal dimensions and nominal volume per chamber.

No. of chambers:	2
Max. no. of shelves:	8
Int. Dim. (WxDxH):	495 x 540 x 500 mm
Ext. dim. (WxDxH):	560 x 580 x 1,150 mm
Max. load shelves:	30 kg
Weight:	36 kg
Temperature resistance:	-20 ... 70 °C
Material:	PMMA
Power supply:	230 V, 50 Hz

Type	Nominal capacity l	Int. width mm	Depth int. mm	Height int. mm	Ext. width mm	Ext. depth mm	Ext. height mm	Min. temp. resistance °C
Control-Maxi 2	156	495	540	500	560	580	1150	-20

Alternativen:
Desiccator Control-Star / Bohlender
Desiccator Control-Big-Star / Bohlender