

Digital refractometer RFM300-M

-High-resolution 4" display with robust push-button keypad

Resolution: 0.001 °Brix0.000001 RI
Accuracy: ±0.01 °Brix (0 ... 30 °Brix)±0.00002 RI (1.32 ... 1.38 RI)
Temperature: 0 ... 70 °C
Power supply: 100 ... 240 V, 50/60 Hz

Type	Measured range	Max. temperature °C	Type of auxiliary energy	PK	Cat. No.
RFM340-M	0 ... 100 °Brix1.32 ... 1.58 RI	70	Mains connection	1	6.322 721

Zubehör:
Flowcell for refractometer series RFM300 / Bellingham + Stanley Ltd.
Funnel flowcell for refractometer series RFM300 / Bellingham + Stanley Ltd.
Printer CMB-910 / Bellingham + Stanley Ltd.
Accessories for printer CMB-910 / Bellingham + Stanley Ltd.
USB thermal printer / Bellingham + Stanley Ltd.
Barcode reader / Bellingham + Stanley Ltd.
USB keyboard / Bellingham + Stanley Ltd.