

Shaking incubator ES-20/80C with cooling

- With Peltier modules for cooling of -12.5 °C below ambient temperature
- Excellent sample temperature uniformity

Please order platforms separately.

Temperature:	RT-12,5 ... 80 °C
Temperature stability:	±0,1 °C
Speed:	50 ... 400 1/min
Shaking stroke:	20.0 mm
Max. loading capacity:	10.6 kg
Dimensions (WxDxH):	620 x 620 x 510 mm
Weight:	50.0 kg
Power supply:	230 V, 50/60 Hz

Type	Min. temperature °C	Max. temperature °C	Min. speed 1/min	Max. speed 1/min	Width mm	Depth mm	Height mm	With timer
ES-20/80C	-12,5	80	50	400	620	620	510	Yes

Zubehör:
Platforms for flasks for shaking incubators series ES-20/80 / Grant Instruments Ltd.
Accessories for shaking incubators series ES-20/80 / Grant Instruments Ltd.

Alternativen:
Shaking incubator ES-20/80 / Grant Instruments Ltd.