

Roller mixer CTR-6

This roller mixer is specifically designed for use in harsh environments such as high CO₂ and humidity and provides reproducible results for cell culture growth. Removable rollers makes the unit flexible and allows for performing various procedures and techniques with various cultivation vessels.

- With remote control to protect the electronics from the CO₂ environment of the incubator
- With stepper motor, guaranteed service life up to 10000 h
- Timer from 1 min to 96 h, continuous
- Stackable with optional stacking kits (up to 3 devices)

Inclination angle:	4 °
Speed:	5 ... 80 1/min
Max. load:	3 kg
No. of rollers:	6
Dimensions (WxDxH):	310 x 262 x 80
With continuous operation:	Yes
Weight:	3 kg
Power supply:	100 ... 240 V·AC, 50/60 Hz

Type	Type of movement	Max. inclination angle °	Min. speed 1/min	Max. speed 1/min	Depth mm	Height	Width mm	Display	PK	Cat.
CTR-6	Rotating/rocking	4	5	80	262	80	310	Digital	1	6.52

Zubehör:
Stacking kits for roller mixer CTR-6 / Grant Instruments Ltd.