

Roller bottles WHEATON^{®}, borosilicate glass, with phenolic resin screw cap

For adherent laboratory cell culture. Meets the requirements of USP Type I and ASTM E-438 Type I, Class A. Autoclavable. Black phenolic resin screw cap with styrene-butadiene insert.

Boro 3.3							
Nominal capacity ml	Ø mm	Height mm	Dimensions (ØxH)	Cap size	With lid	PK	Cat. No.
1380	110	240	110 x 240 mm	38-415	Yes	4	4.662 203
1760	110	285	110 x 285 mm	38-415	Yes	4	4.662 204
3450	110	480	110 x 480 mm	38-415	Yes	4	4.662 205
4230	110	550	110 x 550 mm	38-415	Yes	4	4.662 206
1760	110	270	110 x 270 mm	51-400	Yes	4	4.662 209
2500	110	355	110 x 355 mm	51-400	Yes	4	4.662 210
4230	110	535	110 x 535 mm	51-400	Yes	4	4.662 211

Alternativen:
Roller bottles WHEATON[®], borosilicate glass, with screw cap, PP / DWK Life Sciences Inc.(Wheaton