

Roller mixer

Digital roller mixer with 6 adjustable rollers for tubes and bottles of different sizes. The anti-corrosive, wear-resistant construction enables use in humid and CO₂ environments. A remote control with timer function and a stacking kit are optionally available. Timer 0 to 999 min.

Scope of supply: roller mixer, cord set with various plugs, O-rings

Speed: 1 ... 80 1/min

Max. load: 6.5 kg

No. of rollers: 6

Dimensions (WxDxH): 360 x 378 x 72 mm

Rollers (ØxL): 28 x 259 mm

Permitted ambient temperature: 4 ... 60 °C

With continuous operation: Yes

Weight: 6.5 kg

Rated capacity: 10 W

Plug type: EU/UK/CN

Power supply: 100 ... 240 V, 50/60 Hz

Warranty: 2 years

Type of movement	Min. speed 1/min	Max. speed 1/min	Depth mm	Rollers Ø mm	Height mm	Length rollers mm	Min. permitted ambient temperature °C
Rotating	1	80	378	28	72	259	4