

Weighing boats, PS, square shape, antistatic

White. Light and flat design. Ideal for weighing. Wall thickness approx. 0.2 mm.
Colour: White
Material: PS
Shape: Square

Nominal capacity ml	Width mm	Depth mm	Height mm	PK	Cat. No.
7	45	45	7	500	7.020 016
100	84	84	24	500	7.020 017
250	140	140	21	500	7.020 018

Alternativen:
Weighing boats, PS, square shape, sterile / BRAND
Weighing boats, PS, diamond shape, antistatic / BRAND