

Dual Syringe Continuous Dispenser with Premium Controller Microlab®; 700 Series

Max. capacity: 50 ml
Max. flow rate: 6000 µl/s
Accuracy: ±1 %
Dimensions (WxDxH): 177.8 x 139.7 x 266.7 mm
Weight: 5.9 kg
Power supply: 100 ... 240 V, 50/60 Hz

Type	No. of syringes	Min. capacity µl	Min. flow rate µl/s	Width mm	Depth mm	Height mm	Type auxiliary Main
ML725-CNT Dual Continuous Disp Premium	2	1	0.003	177.8	139.7	266.7	

Zubehör:
Accessories for Microlab® 700 Series / Hamilton Central Europe SRL

Alternativen:
Diluter with Disposable tip hand probe Microlab® 700 Series / Hamilton Central Europe SRL
Dual Syringe Diluter with Advanced controller Microlab® 700 Series / Hamilton Central Europe SRL
Dual Syringe Diluter with Premium controller Microlab® 700 Series / Hamilton Central Europe SRL