

Desiccators Maxi 2, PMMA

Aluminium frame with acrylic panels. Two chambers, two doors, each chamber usable with up to 8 shelves. Including four shelves, two trays and two hygrometers. Internal dimensions and nominal volume per chamber.

- Door with magnetic catch and circular rubber seal
- Four casters (two of the casters with brakes)
- Easy to read electronic hygrometer
- Variable height shelves made of stainless steel
- Telescopic shelves
- Desiccant tray

No. of chambers: 2
Max. no. of shelves: 8
Int. Dim. (WxDxH): 495 x 540 x 500 mm
Ext. dim. (WxDxH): 560 x 580 x 1,150 mm
Max. load shelves: 80 kg
Weight: 34 kg
Temperature resistance: -20 ... 70 °C

Type	Nominal capacity l	Int. width mm	Depth int. mm	Height int. mm	Ext. width mm	Ext. depth mm	Ext. height mm	Min. temp. resistance °C
Maxi 2	156	495	540	500	560	580	1150	-20